




Steril Cerrahi Setler

Sterile Surgical Packs

Disposable Steril El Cerrahi Seti

Disposable Sterile Hand Surgery Pack

Ref. Code
905-001-S



Koli Adet / Unit Box: 8

No/No	Ürün Adı/Product Name	Ebat/Size	Adet/Unit
1	Havlu / Towel	40x40	6
2	Korumalı Cerrahi Önlük / Reinforced Surgical Gown	XL	3
3	Yan Bantlı Örtü / Side Banded Drape	150x240	2
3	Yan Bantlı Örtü / Side Banded Drape	75x90	2
4	Yapışkan Bant / OP Tape	10x50	1
5	Mayo Masa Örtüsü / Mayo Stand Cover	78x145	1
6	Alet Masa Örtüsü / Back Table Cover	150x200	2
7	Sterilizasyon Bohçası / Sterilization Wrap	100x100	1



TSE EN 13795-1
EN ISO 9001:2000
EN ISO 13485:2003






Steril Cerrahi Setler

Sterile Surgical Packs

Disposable Steril El Cerrahi Seti 2

Disposable Sterile Hand Surgery Pack 2

Ref. Code
905-002-S



Koli Adet / Unit Box: 15

No/No	Ürün Adı/Product Name	Ebat/Size	Adet/Unit
1	Havlu / Towel	40x40	6
2	Yan Bantlı Örtü / Side Banded Drape	150x240	2
3	Yan Bantlı Örtü / Side Banded Drape	75x90	2
4	Yapışkan Bant / OP Tape	10x50	1
5	Mayo Masa Örtüsü / Mayo Stand Cover	78x145	1
6	Alet Masa Örtüsü / Back Table Cover	150x200	2
7	Sterilizasyon Bohçası / Sterilization Wrap	100x100	1



TSE EN 13795-1
EN ISO 9001:2000
EN ISO 13485:2003

