




## Steril Cerrahi Setler

### Sterile Surgical Packs

#### Disposable Steril Plastik Cerrahi Seti

#### Disposable Sterile Plastic Surgery Pack

Ref. Code  
915-001-S



Koli Adet / Unit Box: 20

No/No	Ürün Adı/Product Name	Ebat/Size	Adet/Unit
1	Yan Bantlı Örtü / Side Banded Drape	80x90	2
2	Bantlı Ayak Örtüsü / Foot Drape	150x180	1
3	Anestezi Örtüsü / Anesthesia Drape	150x240	1
4	Alet Masa Örtüsü / Back Table Cover	150x200	1
5	Sterilizasyon Boşçası / Sterilization Wrap	50x50	1



TSE EN 13795-1  
EN ISO 9001:2000  
EN ISO 13485:2003

