



Steril Cerrahi Setler

Sterile Surgical Packs

Katater Set

Vesical Catheterization Pack

Ref. Code
919-001-S



Koli Adet / Unit Box: 40

No/No	Ürün Adı/Product Name	Ebat/Size	Adet/Unit
1	Düz Örtü DL / Plain Drape	35x40	1
2	Delikli Örtü DL (4.5x4.5 dlk.) / Drape With Aperture	40x45	1
3	Gaz Bez Spanç / Gauze Sponge	10x10	3
4	Gaz Top / Gauze Ball	st	3
5	Enjektör / Syringe	10 cc	1
6	Forseps / Forceps	st	1
7	Koter Forseps / Cautery Forceps	st	1
8	Eldiven / Gloves	m	1
9	Jel / Gel	20 g	1
10	Genital Mendil / Genital Wipe	st	1
11	Tepsi / Tray	st	1



TSE EN 13795-1
EN ISO 9001:2000
EN ISO 13485:2003

